

## Universal Neutralizing Sorbent #UNS4 - 3lb Shaker, 4 each/package



Neutralizes, solidifies and indicates safe pH range in one easy step.

### UNS4 Technical Information

#### Warnings & Restrictions:

##### UNS Notice

**SAFE TO USE!** Reactions with aqueous solutions and acids will generate some heat depending on the strength of the solution. Reactions with bases generate little or no heat. As a general rule of thumb, apply the same amount of weight as the spill (1 lb. to 1 lb.). Example: 98% sulfuric acid is 14.97 lbs. per gallon, apply 15 lbs. of UNS to 1 to 1 ratio. This rule of thumb varies based on the specific chemical and its strength. **CAUTION:** Only qualified, trained personnel should attempt the cleanup of a hazardous spill. Wear appropriate personal protective equipment (PPE). Product may cause irritation to the eyes and skin. Some heat may be generated when used on strong acids or bases. For eye contact, flush eyes with water, seek medical attention if necessary. For skin contact, wash with soap and water. For inhalation of vapor, move to fresh air and seek medical attention if necessary. Disposal must be in compliance with local, state and federal regulations.

#### Technical Documents:

(Available at [spilltech.com](http://spilltech.com))

Product Data Sheet (PDS)

Safety Data Sheet (SDS)