

Multi-Purpose Neutralizing Spill Kit #MPN-KIT - 12" dia. x 10", 1 each/package



Treat a wide range of hazardous chemicals and solvents.

- Safely clean spills of sulfuric, hydrochloric, nitric and phosphoric acid, hydrofluoric acid, caustics and common hydrocarbon solvents all with one kit
- Granular formulas require less clean-up and do not create "dust clouds"
- Kit contents include: (2) 1.5 lb shakers of Ultimate Acid Eater Neutralizing Absorber, which contains, neutralizes and solidifies a variety of acids, including sulfuric, hydrochloric, nitric, or phosphoric; (1) 1.5 lb shaker of HF Acid Eater Neutralizing Absorbent, which contains, neutralizes and solidifies hydrofluoric acid; (1) 1.5 lb shaker of Base Eater Neutralizing Absorbent, which contains, neutralizes and solidifies caustics; (1) 1 lb shaker of BioRem-2000 Solvent Absorber, which absorbs spills of many common hydrocarbon solvents; 2 Goggles, 2 Aprons, 2 Pairs of Gloves; Disposal Bag & Tie; Broom & Dustpan set
- Kit packaged in a 3.5 gal portable pail

MPN-KIT Specifications

Dimensions:	ext. dia. 12" x 10" H
Shipping	12" W x 12" L x 10" H
Dimensions:	
Sold as:	1 per package
Composition:	Various
Weight:	15 lbs.
# per Pallet:	48
Incinerable:	No
Ship Class:	55

Metric Equivalent Specifications

Dimensions:	ext. dia. 30.5cm x 25.4cm H
Shipping	30.5cm W x 30.5cm L x 25.4cm H
Dimensions:	
Weight:	6.8 kg

MPN-KIT Technical Information

Warnings & Restrictions:

There are no known warnings and restrictions for this product.

Technical Documents:

(Available at spilltech.com)

Product Data Sheet (PDS)

Safety Data Sheet (SDS)



SpillTech

Online:
spilltech.com

Phone:
1-800-228-3877 (N.Am.)
1-770-475-3877 (Other)

Fax:
1-800-872-3764 (N.Am.)
1-770-410-1812 (Other)

Email:
sales@spilltech.com